

# SMART REAGENTS

## Magnesium Chloride Reagent (MgCl)



### Features

- Formulated from analytical grade reagent
- Consistency from lot to lot
- Dissolve and go for greater convenience
- Eliminates weighing errors

### Product description

Magnesium chloride is an important coagulant and a powerful stimulant of the immunological system. The compound is used in food as a nutritional supplement and in industrial processing. Medicago's MgCl is more than 99% pure and supplied as pre-weighed powder in sealed pouches, giving 1000 ml of 1 M magnesium chloride.

### Applications

- Coagulant
- Immune stimulant
- Nutritional supplement
- Industrial processing

### Directions for use

Empty one pouch of magnesium chloride reagent in a laboratory flask or beaker placed on a magnetic stirrer. Add 300 ml of deionized water and stir for a few minutes. Adjust the volume up to 1000 ml, stir until full dissolution and the solution is ready to use.

### Shipping and storage

The magnesium chloride is shipped at room temperature. Store the pouches in a dry place at room temperature. Shelf life is three years.

### Specifications

Chemicals	Analytical grade
Format	Exactly pre-weighed powder
Volume	1000 ml
Shelf life	Three years after production date

### Tips and hints

If the contents of the pouch is not properly dissolved, make sure:

- the water temperature is 25°C (do not exceed this temperature)
- the solution is properly stirred.

### Certifications

Medicago's laboratories and manufacturing site in Uppsala are ISO 9001:2008 and ISO 13485:2003 certified. Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



### Ordering information

Article no.	Product name	Pack size	Solution vol.
12-9180-5	Magnesium Chloride, 1 M	5 pouches	1000 ml/pouch

v. 01 2010-06