



Immunosorb A

Product description

Immunosorb A is a protein A immobilized on cross-linked agarose matrix. For purification of human, mouse and rabbit immunoglobulins. Contains 20 % ethanol as preservatives. For laboratory use only.

Form

Immunosorb A is supplied as a liquid substrate.

Instructions for use

Immunosorb A is a ready to use substrate. There is no additional mixing or reagents required. It is not necessary to warm the product to room temperature before use.

Product specifications

APPEARANCE	White colour
PROTEIN CONTENT	Protein A is recombinant and is derived from the surface of a bacterium and consists of non-hazardous material. It has been cloned and expressed in E. coli and is >95% pure.
STRUCTURE	The protein binds to specific classes of immunoglobulins.
NATURE	The material is non-cultural and non-infectious and is not pathogenic for animals or plants.
OTHER INFORMATION	The material does not contain any animal or cell culture derived products or additions such as albumin or serum. The material does not contain any toxic or harmful substances. The substances are non-toxic and non-hazardous and will be used for research purpose only.

Ordering information

PRODUCT NUMBER	PRODUCT DESCRIPTION	SOLUTION VOLUME
10-1257-200	Immunosorb A	200 ml
10-1257-100	Immunosorb A	100 ml
10-1257-50	Immunosorb A	50 ml

Revision 01