



IMMUNO REAGENTS

Carbon suspension for CIA



Features

- Ready-to-use
- Specially developed carbon suspension for CIA

Product description

Medicago's Carbon Suspension is specially developed for Carbon Immunoassays (CIA). The suspension contains microscopic carbon particles that attach to IgG molecules in many different mammalian sera, including human serum.

For laboratory use only.

Applications

- Custom carbon Immunoassay studies
- Detection of *E. cuniculi* and *T. gondii* in CIA

Shipping and storage

The product is shipped at room temperature and is to be stored at 2°C to 8°C.

Certifications

Medicago's laboratories and manufacturing site in Uppsala are ISO 9001:2008 and ISO 13485:2003 certified. Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



Specifications

Content	Carbon particles in suspension
Appearance	Liquid solution
Storage	2°C to 8°C

Ordering information

Article no.	Product name	Pack size
18-6010-3	Carbon suspension	3 ml

Directions for use

The product is a ready-to-use suspension. Use lab coat and gloves. If spill or leak occurs, clean surface immediately.

v. 01 2010-08